

# Convenient COVID Ag Testing at the Point of Need



CovIDx™

## Introducing the CovIDx™ SARS CoV-2 Rapid Antigen Test

FOR PROFESSIONAL USE ONLY\*




**RAPID RESULTS:** facilitate immediate patient isolation and prioritization of confirmatory testing

**FLEXIBLE WORKFLOW:** pre-filled extraction vials and individually wrapped reagents enable multiple independent workstations for increased test distribution

**EASY-TO-USE:** instrument-free, user-friendly test procedure for non-lab settings


**CONVENIENT:** all materials included in the kit for hassle-free logistics

**ACCURATE:** 100% PPA with PCR for samples with Ct values  $\leq 30^1$




Point-of-Care  
TEST RESULTS  
within  
**20**  
MINUTES

### SIMPLE PROCEDURE




### SIMPLE RESULT INTERPRETATION




### SARS-CoV-2 Viral Load Over Course of Infection<sup>2</sup>

Frequent testing with antigen tests can identify people when their infection is most likely to be transmissible.<sup>3</sup>



### CoviDx Detection Compared to PCR Ct Values<sup>1</sup>

Increased cycle threshold (Ct) indicates decreasing sample viral concentration.



### CLINICAL PERFORMANCE

#### Prospective Data (Nasal Swabs)<sup>1</sup>

CoviDx SARS-CoV-2 Rapid Antigen Test	PCR test					
	Ct ≤ 30			All		
	Positive	Negative	Total	Positive	Negative	Total
Positive	31	5	36	41	5	46
Negative	0	96	96	9	96	105
<b>Total</b>	<b>31</b>	<b>101</b>	<b>132</b>	<b>50</b>	<b>101</b>	<b>151</b>
Positive Percent Agreement (PPA) Sensitivity	<b>100.0%</b> (95% CI: 89.0% - 100.0%)			<b>82.0%</b> (95% CI: 69.2% - 90.2%)		
Negative Percent Agreement (NPA)	<b>95.0%</b> (95% CI: 88.9% - 97.9%)			<b>95.0%</b> (95% CI: 88.9% - 97.9%)		

#### Retrospective Data (Nasopharyngeal Swabs)<sup>1</sup>

CoviDx SARS-CoV-2 Rapid Antigen Test	PCR test					
	Ct ≤ 30			All		
	Positive	Negative	Total	Positive	Negative	Total
Positive	33	0	33	33	0	33
Negative	0	29	29	4	29	33
<b>Total</b>	<b>33</b>	<b>29</b>	<b>62</b>	<b>37</b>	<b>29</b>	<b>66</b>
Positive Percent Agreement (PPA)	<b>100.0%</b> (95% CI: 89.6% - 100.0%)			<b>89.2%</b> (95% CI: 75.3% - 95.7%)		
Negative Percent Agreement (NPA)	<b>100.0%</b> (95% CI: 83.3% - 100.0%)			<b>100.0%</b> (95% CI: 83.3% - 100.0%)		

1. CoviDx [package insert] PM-128.0. Sarasota, FL: Lumos Diagnostics; 2021. 2. Adapted from Crozier A. Put to the test: Use of rapid testing technologies for Covid-19. BMJ. 2021;372:n208. 3. Mina MJ, Parker R, Larremore DB. Rethinking Covid-19 test sensitivity - A strategy for containment. N Engl J Med. 2020;383:e120.

#### For use under Health Canada Interim Order Authorization only

- Refer to package insert for full performance, warnings, and limitations
- PCR is the gold standard for diagnosing COVID-19
- Testing should follow provincial or territorial guidance and be administered regularly to employees in accordance with the interim guidance on the use of rapid antigen detection tests for identifying SARS-CoV-2 infection