

# Resona I9T Diagnostic Ultrasound System

## Performance Specifications

### System Overview

#### Advantages of ZST<sup>+</sup> platform

- Advanced Acoustic Acquisition
- Dynamic Pixel Focusing (DPF)
- Sound Speed Compensation (SSC)
- Total Recall Imaging (TRI)
- Powerful Processing Architecture
- Enhanced Channel Data Processing
- ZONE Sonography<sup>®</sup> Technology
- Up to 8,257,536 channels

#### Application

- Abdomen
- Obstetrics
- Gynecology
- Cardiology
- Small parts
- Urology
- Vascular
- Pediatrics
- Nerve
- Emergency & Critical
- Pelvic Floor
- Others

#### Transducer Types

- Curved array transducer
- Linear array transducer
- Phased array transducer
- 4D Volume transducer
- Pencil transducer

#### Advanced Imaging Techniques

- THI (Tissue Harmonic Imaging) and PSH<sup>™</sup> (Phase Shift Harmonic Imaging)
- iBeam<sup>™</sup> (Spatial Compound Imaging)
- iClear<sup>™</sup> (Speckle Suppression Imaging)
- iClear+
- iTouch<sup>™</sup> (Auto Image Optimization)
- Echo Boost<sup>™</sup>
- Zoom/iZoom (Full Screen Zoom)
- FCI (Frequency Compound Imaging)
- B steer
- ExFOV (Extended Field of View)
- HD Scope
- Smart Track
- HR Flow<sup>™</sup> (High Resolution Flow)
- Glazing Flow

#### Imaging Modes

- B-Mode
- M-Mode/Color M-mode
- Color Doppler Imaging
- Power Doppler
- Imaging/Directional PDI
- Pulsed Wave Doppler
- Continuous Wave Doppler
- Free Xros M<sup>™</sup> (Anatomical Mmode)

Free Xros CM<sup>™</sup> (Curved Anatomical M-mode)

#### Function Modes

- iScape<sup>™</sup> View (Panoramic Imaging)
- TDI (Include TVI, TVD, TVM, TEI)
- TDI QA (TDI Quantitative Analysis)
- UWN+ (Ultra Wideband Non-linear Plus) Contrast Imaging<sup>™</sup>
- Contrast Imaging QA (Quantitative Analysis)
- LVO (Left Ventricular Opacification)
- Volume CEUS
- Natural Touch Elastography Imaging
- STE Imaging (Sound Touch Elastography)
- STQ Imaging (Sound Touch Quantification)
- Endocavity STE
- High frame rate STE
- Stress Echo
- TTQA (Tissue Tracking Quantitative Analysis)
- Smart 3D<sup>™</sup> (Freehand 3D)
- Real-time 4D
- iPage+ (Multi-Slice Imaging)
- SCV+ (Slice Contrast View)
- Color 3D
- iLive
- DICOM
- Clinical Measurement Package
- Smart Pelvic Floor
- Smart OB<sup>™</sup> (Auto OB measurement)
- Smart Fetal HR (Fetal Heart Rate)
- IVF
- IMT
- Auto EF
- V-Mapping
- iNeedle<sup>™</sup> (Needle Visualization Enhancement)
- ECG function

#### Other Features

- Ultrasound gel
- Ultrasound gel warmer
- Built-in wireless adapter
- Replaceable battery assembly
- Central brake
- Probe adapter
- DVD assembly
- Wipes box bracket
- DVR
- Ambient light
- iStorage
- iWorks<sup>™</sup> (Auto Workflow Protocol)
- MedSight
- MedTouch
- UltraAssist (Off-line software)
- UltraView<sup>™</sup> (Off-line analysis software)
- Touch gestures
- Anti-virus software: McAfee, Window
- iVocal



#### Language Support

- Software: English
- Keyboard input: English

#### Physical Specifications

##### Dimension and Weight

The control panel and the monitor are in the lowest position  
Configured with floating support arm and 23.8 inch monitor

- Depth: 1020±20 mm;
- Width: 550±10 mm;
- Height: 1000±20 mm
- Weight: 105 kg±4 kg (net weight, standard configuration but not including the probe)

#### Audio Speaker

Stereo audio speakers

#### Dual-Wing Floating Support Arm

- Rotate angle: 90±5 degrees (to the left); 150±5 degrees (to the right)
- Tilt angle (when positioned vertically): 20±5 degrees (backward); 85±5 degrees (forward)
- From front to back: 300±20 mm
- From bottom to top: 150±20 mm

#### Wheels

- Diameter: 125 mm
- When the central brake is configured: Central brake for total lock and break

## Performance Specifications

### Physical Specifications (continued)

When the central brake is not configured: 3 castors for total lock and break, and 1 castor for direction lock and break

### Transducer Port and Holder

Transducer ports with dust prevention: 5 active ports and 1 pencil probe port

Support hot plug with active indicator lights

Transducer holder: 5, plus 1 dedicated endocavity transducer holder and 1 dedicated pencil transducer holder

### Electrical Power

Voltage: 100-240 V~

Frequency: 50/60 Hz

Power consumption: 650 VA

### Operating Environment

Ambient temperature: 0-40° C

Relative humidity: 20%-85% (no condensation)

Atmospheric pressure: 700 hPa-1060 hPa

### Storage & Transportation Environment

Ambient temperature: -20-55° C

Relative humidity: 20%-95% (no condensation)

Atmospheric pressure: 700 hPa-1060 hPa

### System Noise

≤26 dB @25° C

### User Interface

#### Floating Control Panel

Brightness adjustable for the backlight of the whole control panel

Full-sized, backlit QWERTY keyboard

iConsole: intelligent control panel for clinical-exam specific layout and adaptive adjustment: 6 programmable E-ink keys for dynamic display of user-defined functions

Automatic light indication for residual battery power

Full-space floating control panel adjustment and can be fixed at any position (when centered in the trackball):

Left/right rotation: 180±5 degrees (90 degree for both left and right)

Left/right adjustment: 1100 mm±50 mm (550 mm for both left and right)

Down/up adjustment: 300±20 mm

Front/back adjustment: 350±20 mm

### Monitor

23.8-inch bezel-less LED monitor with high resolution

Resolution: 1920x1080

Viewing angle: 178 degrees

Digital on screen display of brightness and contrast controls

Automatic adjustment of monitor light with the changing environment

Automatic LED brightness

Tilt/Rotate independent adjustment

Tilt angle range: 105 degrees

Rotate angle range: 240 degrees

### Touchscreen

15.6-inch high sensitivity antiglare color touch screen

Resolution: 1920\*1080

Digital brightness and contrast adjustment through preset

Angle adjustable range: 40 degrees

Viewing angle: 170 degrees

Support touch screen gestures

Support either hand writing or with gloves on

Editable buttons: long press to add, delete or move buttons

Moveable 3D/4D tabs

Clinical scenario-based 3D/4D user interface

Digital TGC

Short-cut switch of latest used probe & exams

### Touch Gestures

Swipe down/up: display/remove projected image on touch screen

Swipe horizontally: page up/down or review images/cine loops one by one

Swipe from left edge to right: display hidden menu on projected image

Image parameter adjustment

Measurement on projected image on touch screen

Zoom in/out the projected image on touch screen

Rotate or erase on projected 3D/4D image on touch screen

8 user defined gestures using two fingers for more functions, such as freeze, print, activate specific imaging modes, measurements, and some other special functions

### System Boot-Up

Boot-up from shut-down: <30 sec

Boot-up from stand-by: <5 sec

Shut-down: <30 sec

### Comments

Supports text input and arrow

Voice annotation: record voice as annotation for images and cine

Support freehand marking on touch screen

Adjustable text size and arrow size

Supports home position

Covers various application

User customizable

### Bodymark

More than 232 bodymarks for versatile application

User customizable

### Imaging Parameters

#### Overview

Echo-enriched Beamforming

Up to 55296 channels

Up to 8-beamforming

#### B-Mode

Display formats:

Single (B)

Dual (B+B), support by B/ M/ Color/ Power/ PW CW/ Color M mode

Quad (4B), support by B/ Color/ Power

iClear™/ iClear+™: Off; 1-7, 7 steps

iBeam™: Off, 3 steps

iTouch™: On/off

FCI (Frequency Compound Imaging)

Dual Live: On/off

Imaging Parameters

#### Imaging

Image quality: Pen/Gen/Res/HPen/HGen/HRes/HGen-FFR/HRes-FFR (dependent on transducer)

B steer: 5 levels, available on linear transducers

ExFOV: off, 1-2 (dependent on transducer)

Depth: 30 levels, 1.5-40 cm

Frame rate (max): 1701 f/s

Acoustic output power: dependent on transducer

TGC/LGC: 8 segments on touch screen

Dynamic range: 30-260 (dependent on transducer)

Gain: 0-100, 1/step

FOV: continuously adjustable

Line density: L/M/H/UH

Persistence: 0-7 levels

Horizontal Scale: on/off

L/R flip and U/D flip: on/off

Rotation: 0°, 90°, 180°, 270°

TSI: general/muscle/fluid/fat

Gray Map: 8 types

Tint map: off, 8 types

Echo Boost: On/Off

Smooth: 0-6 levels

HD Scope: off, 1-3 levels

SSC (Sound Speed Compensation): On/Off

Free view: -45°~45°, 5°/step

Dehaze: 0-30 levels

Ref Lines: on/off (under GYN and Pelvic Floor exam mode and using intracavity probe)

## Performance Specifications

### Imaging (continued)

V 1:1:	on/off (available on linear probe under dual-split mode)
ExtImage:	On/Off
Auto Merge:	On/Off
ZoneVue	
Edge enhance:	0-6, 1/step

### THI and PSH™

Patent PSH™ technology, obtains purer harmonic, better contrast resolution, higher SNR, exceptional high frequency harmonic

iClear™ available

Image quality: HPen/HGen/HRes or HPen/HPen-FFR/HGen/HRes/HRes-FFR (depends on transducers)

Echo Boost™: On/Off

### M-Mode

Display formats: V2:3, V3:2, V3:1, H2:3, Full (V: vertical, H: horizontal)

Color M-mode available

Acoustic output power:	same as B
Depth:	same as B
Dynamic range:	30-180, 5/step
Gain:	0-100, 1/step
M sweep speed:	6 steps
M soften:	0-4, 1/ step
Tint map:	off, 8 types
Gray Map:	8 types
Edge enhancement:	0-3, 1/ step

### Color Doppler Imaging

Dual live:	On/Off
HR Flow™:	High Resolution Flow provides better image quality and flow sensitivity
Image quality:	Pen/Gen/Res (color), 1 level (HR Flow)
Max velocity:	146.5 cm/s
Steer:	available on linear transducer
Max frame rate:	529 f/s
Acoustic output power:	same as B mode
Gain:	0-100, 2/step
ROI size/position:	adjustable
Scale:	max. 30 steps
Baseline:	-8 – 8, 1/step
Wall filter:	8 steps, 5-433 Hz
PRF:	0.1-12.6 kHz
Packet size:	0-3, 1/step
Flow state:	L/M/H
Smooth:	0-6, 1/step
B/C align:	On/Off
Priority:	0%-100%, 1%/step
Color map:	V0-V10; VV0-VV9
Invert:	On/Off
Auto Invert:	On/Off
Persistence:	0-6, 1/step
Velocity tag:	On/Off
Line density:	L/M/H/UH

iTouch™:	on/off
Smart track:	On/Off
Glazing Flow:	On/Off, L/M/H

### Power Doppler Imaging

Dual live:	On/Off
HR Flow™:	High Resolution Flow provides better image quality and sensitivity Support directional power Doppler
Image quality:	Pen/Gen/Res (Power), 1 level (HR Flow)

Acoustic output power:	same as B
Dynamic range:	10-70, 5/step
Gain:	0-100, 2/step
ROI size/position:	adjustable
Steer:	available on linear transducers
Scale:	max. 30 steps
Wall filter:	8 steps
PRF:	max. 12.6 kHz
Packet size:	0-3, 1/step
Flow state:	L/M/H
Smooth:	0-6, 1/ step
B/C align:	On/Off
Priority:	0%-100%, 1%/step
Color map:	4 types
Directional color map:	4 types
Persistence:	0-6, 1/step
Line density:	L/M/H/UH
Invert:	On/Off
iTouch™:	On/Off
Smart track:	On/Off
Glazing Flow:	On/Off, L/M/H

### PW/CW-Mode

Display formats:	V2:3, V3:2, V3:1, H2:3, Full, Duplex/Triplex (PW only) (V: vertical, H: horizontal)
Image quality:	Pen/Gen/Res
PW velocity:	max. 8.681 m/s min. 2.03 cm/s
CW velocity:	max. 39 m/s min. 5.34 cm/s
Sample volume size:	0.5-30 mm (PW only)
Sample gate depth:	continuously adjustable
Baseline:	9 steps
PW Steer:	available on linear transducer
Volume:	0%-100%, 2%/step
PW PRF:	1.0-23.1 kHz
CW PRF:	0.3-104.0 kHz
Gain:	0-100, 2/step
Dynamic range:	24-72, 2/step
Sweep speed:	6 steps
Wall Filter:	0 ~ 9 steps
Invert:	On/Off
Auto invert:	On/Off
Angle correction:	-89~89 degrees, 1/step
Quick angle:	0, -60, 60 degrees
Gray map:	10 types
Tint map:	Off; 8 types

Time/frequency resolution:	0-6, 1/step
HPRF:	On/Off
Auto calc:	On/Off
Auto calc cycle:	1, 2, 3, 4, 5
Auto Calc Loop	
Trace Sensitivity:	0-5, 1/step
Trace Smooth:	Off, 1-4, 1/step
Trace area:	Above, Below, All

### Free Xros M™

Display formats: V2:3, V3:2, V3:1, H2:3 (V: vertical, H: horizontal)

Color Free Xros M available

Up to 3 lines	
Display all lines	
Sweep speed:	6 steps
M Tint map:	Off; 8 types
Gray Map:	8 types

### Free Xros CM™

Only available in TDI mode

Display formats: V2:3, V3:2, V3:1, H2:3 (V: vertical, H: horizontal)

Sweep speed:	6 steps
Tint map:	Off; 8 types
Gray Map:	8 types
Edit, Undo, Delete function for curved line	

### iBeam™

Spatial compound imaging  
Off, 1-3, 1/step

### iClear™

Speckle suppression imaging  
Available on B, 3D, 4D mode

### iTouch™

Auto image optimization	
B-mode:	gain, TGC, dehaze
Color:	gain, color box position
Power:	gain
PW:	gain, scale, PRF, WF, SV size, SV position, steering angle
Contrast imaging:	gain

### Echo Boost™

Available in cardiac exam mode when using a phased array probe  
Improve the homogeneity through the whole field of view  
Better noise control in cardiac chambers and muscles

### Zoom

Zoom:	Spot zoom (write zoom), Pan zoom (read zoom) 0.8x-10x
iZoom:	convertible 3 steps; normal image, zoom standard area, zoom only image area

### QSave

Quickly save image parameter setting after image adjustment done  
Support Save, Create, Restore

## Performance Specifications

### Imaging (continued)

#### IP (Image Process) Preset Manager

Support save, create, restore or delete IPs  
Quickly switch to different exam modes without returning to B mode  
A default set of image parameters under each IP are provided for different exam modes  
Only image parameters are switched while measurements, comments, and bodymarks remain the same

#### 3D/4D Preset Manager

The scenarios and subpresets can be renamed, restored, deleted, added, set to active, or moved  
Show scenario and subpreset parameters  
Provide multiple groups of preset 3D/4D parameters based on different application scenarios to quickly obtain expected image effect

### Tissue Velocity/Energy Imaging (included in TDI option)

Available on phased array, and SC6-1s, SC8-2s, SP5-1s, and P7-3Ts  
Dual live: side by side displays B and B+TVI  
Max frame rate: 3175 f/s  
PRF: 0.4-14.9 kHz

Acoustic output power: same as B  
Gain: 0-100, 2/step  
Dynamic range: 10-70, 5/step (TEI only)  
ROI size/position: adjustable  
Scale: max. 30 steps  
Wall filter: 8 steps  
Packet size: 0-3, 1/step  
Flow state: L/M/H  
Smooth: 0-6, 1/step  
B/C align: On/Off  
Priority: 0%-100%, 1%/step  
Color map: 10 types (TVI), 8 types (TEI)  
Invert: On/Off  
Persistence: 0-6, 1/step  
Velocity tag: On/Off (TVI only)  
Line density: L/M/H/UH  
Image quality: 2 levels

### Tissue Velocity Doppler (included in TDI option)

Display formats: V2:3, V3:2, V3:1, H2:3, Full, Duplex/Triplex (V: vertical, H: horizontal)  
Sample volume size: same as PW  
Sample gate depth: adjustable  
Sample volume depth: continuously adjustable  
Scale: 30 levels  
Volume: 0%-100%, 2%/step  
PRF: 0.7-23.1 kHz  
Gain: 0-100, 2/step  
Baseline: 9 steps  
Dynamic range: 24-72, 2/step  
Sweep speed: 6 steps

Wall Filter: 10 steps  
Invert: On/Off  
Angle correction: -89-89 degrees, 1/step  
Quick angle: 0, -60, 60 degrees  
Gray map: 10 types  
Tint map: Off; 8 types  
Image quality: 2 levels  
Time/frequency resolution: 0-6, 1/step

### Tissue Velocity Motion (included in TDI option)

Display formats: V2:3, V3:2, V 3:1, H2:3, Full (V: vertical, H: horizontal)  
Acoustic output power: same as B  
Gain: 0-100, 2/step  
M sweep speed: 6 steps  
Smooth: 0-6, 1/ step  
Color Map: 10 types  
Image quality: 2 levels  
Persistence: 0-6, 1/ step  
Packet size: 0-3, 1/ step  
Priority: 0%-100%, 1%/step  
Velocity tag: On/Off  
Tissue state: L/M/H

### Smart 3D™

Smart 3D  
Acquisition preparation:  
3D/4D Scenario setting: Routine  
Acquisition Methods: Rocked, Linear  
Reset VOI: On/Off  
Flip VOI: On/Off  
Angle: 10-80°  
Distance: 10-200 mm  
Acquiring Time: 1.0s-20.0s  
VR:  
3D/4D Scenario: Routine (Surf., iLive Gen., iLive Transp., Skeleton)  
Reset: All, Orientation, Curve  
VOI: On/Off/Fixed  
Active quadrant: A, B, C, VR  
VR orientation: 0°, 90°, 180°, 270°  
Flip: flip VR  
Sync: synchronize VR with selected plane  
Orientation Assist: On/Off  
Threshold: 0-100%, 1%/step  
Opacity: 0-100%, 5%/step  
Brightness: 0-100%, 2%/step  
Contrast: 0-100%, 2%/step  
Smooth: 0-10, 1/ step  
Depth VR: Off/Black/Cyan/Blue/ Rose  
Tint: off; 8 types  
Degree: 10-80°  
Distance: 10-200 mm  
MPR:  
Active quadrant: A, B, C  
Gray Map: 1-8  
Brightness: 0-100%, 2%/step

Contrast: 0-100%, 2%/step  
iClear: Off; 7 types  
Tint: Off; 8 types  
Thickness: 0-30 mm  
Adv.:  
Direction: Up/Down, Left/Right, Front/Back, Down/Up, Right/Left, Back/Front  
VR Refine: Off; 7 steps  
Surface enhance: 0-7, 1/step  
MagiClean: Off/Low/Mid/High/Max  
Inversion: On/Off  
A3:1: On/Off  
Move light: On/Off  
Degree: 10-80°  
Distance: 10-200 mm  
Main render: Surface, Max, Min, X Ray, iLive  
Sub render: Surface, Max, Min, X Ray  
Mix: Set the mix ratio of the two render modes

### Tool

Edit:  
Rubber: On/Off  
Eraser Diameter: 8-80, 1/step  
Cut (area selection): Polygon, Contour, Rectangle, Line  
Undo: Undo, Undo All  
3D Layout:  
Niche Views: Inner, Outer  
3 Slice  
Active Quadrant: A, B, C,  
Niche/3 Slice  
Auto rotation:  
Position: Set Start/Set End  
Direction: Left/Right, Up/Down  
Step: 1-15°  
Quick Angle: 30-180°  
Rotation control: play, single loop, loop  
Save AVI to USB

### 4D

Available on all volume transducers  
Static 3D and real time 4D  
Acquisition preparation:  
4D frame rate: max. 80 vps  
3D/4D Scenario setting: Smart Scene3D (Spine, Brain, Long Bone, Face, Endometrium, Pelvic), Routine, iLive Pro, Bone, Tissue (not all scenarios are listed)  
Refresh: On/Off  
Angle: 10-120°  
Quality: low1, low2, mid, high1, high2  
Other parameters are the same as Smart 3D  
VR:  
3D/4D Scenario: Smart Scene 3D, Routine, iLive Pro, Bone, Tissue, Routine (not all scenarios are listed)  
3D iClear: Off; 7 steps  
Face+: Off, 3 steps

## Performance Specifications

### Imaging (continued)

Auto Play:	Stop, x1, x2, x3, x1/2, x1/3
Frame:	Select a frame
Other parameters are the same as Smart 3D	
MPR	
3D iClear:	Off; 7 types
Other parameters are the same as Smart 3D	
Adv.:	The parameters are the same as Smart 3D
Tool	
The parameters are the same as Smart 3D	
3D Reference Point	
Enable the operator to define one or more reference points on MPRs, which are then projected to VR image; helpful for the operator to better understand the corresponding spatial relations of VR image and MPRs	
Display:	Point only, H line, V line
Delete All	
Hide All	
Color 3D	
Available on volume transducers Supports Color and Power mode	
Only Available in Static 3D mode	
iPage+	
Slice display mode:	Slice only, Slice with SCV
RefMPR:	A Plane, B Plane, C Plane
Reset Ori:	On/Off
Sync MPR:	A/B/C
Quick Rotation:	90°, -90°, -180°
Slice layout:	2*2, 3*3, 4*4, 5*5
Slice Number:	3-25
Spacing:	0.5-10mm
Line (Start)/Line (End)	
Single window	
Display Ref plane	
Hide Ref plane	
Vertically display ref plane	
Horizontally display ref plane	
Move plane	
SCV+	
Reset:	All
Orientation Assist:	On/Off
Active Quadrant:	A, B, C
Main render modes:	Surface, Max, Min, X-ray
Sub render modes:	Surface, Max, Min, X-ray
Mix:	Set the mix ratio of the two render modes
CMPR	
Reset Curve	
Undo Last	
Flip:	Thickness, Direction
Ref. Image (A, B, C), 1, 2, 3	
Trace Options:	Line, Trace, Spline

Quick Rotation:	90°, -90°, -180°
Hide Annotation	
Rotate RL	
Support labeled measurements	
iLive	
Shading:	0-10, 1/step
Grad View	
Hyaline:	On/Off, 0~100%, 5%/step
Move Light	
Light Position:	6 default positions selectable
Light 1/2/3:	Parallel, Point, Torch
VL Saturation:	0~100%, 1%/step
VL Hue:	0~100%, 1%/step
VL Distance:	0.0 – 5.0
VL Bright:	0~100%, 1%/step
VL Angle:	0~100%, 1%/step
Tint:	0-4
Reset	
Classic/IntPoint/ExtPoint/Parallel/ Torch/ 3-Light/ User 1/User 2	
Copy to:	Copy the current lighting mode to customized lighting mode "User 1" or "User 2"

### Smart Track

Available on linear probes	
Enable the function under Color/Power mode, the angle and the position of the ROI are adjusted automatically	
Enable the function under Color/Power+PW mode, the angle and the position of the PW sampling line, SV size, SV angle and SV position are adjusted automatically	

### iScale™ View

Acquisition method:	B and Power
Supports speed indicator	
Actual size:	On/Off
Fit size:	On/Off
Ruler:	On/Off
Tint map:	Off; 8 types
Rotation:	0~355 degrees, 5/step

### iNeedle

Needle visualization enhancement	
Available on linear SC6-1s, and C6-2Gs probes	
Needle direction:	Auto, Left, Right
B/iNeedle:	On/Off
In-plane biopsy and Out-plane biopsy	
GPS real-time guiding biopsy	

### UWN+ Contrast Imaging™

UWN+ (Ultra Wideband Non-linear Plus) contrast imaging technology, which provides exceptional contrast agent detecting capability, not only extracts second harmonic, but also nonlinear fundamental signals	
Available on SC6-1s, SC8-2s, C11-3s, C6-2Gs, L9-3s, L14-3Ws, L20-5s, V11-3Hs, SP5-1s, SD8-1s, DE11-3Ws, and ELC13-4s probes	
Micro Flow Enhancement (MFE) available	

Timer1:	On/Off
Timer2:	On/Off
Pro capture:	captures prospective image
less than 480 s preset table	
Retro capture:	captures retrospective image
less than 120 s preset table	
Dual live:	side by side displays tissue
image and contrast image	
MFE:	On/Off
MFE period:	0.1 s, 0.2 s, 0.4 s, 0.6 s, 0.8 s,
1.0 s, MAX	
Destruct:	instantly destroy contrast
bubbles	
Destruct AP:	-30~0 dB, 0.3/step
Destruct time:	500-2000 ms
Contrast Agent:	SonoVue and SonaZoid
iClear:	Off; 7 steps
Mix:	mix contrast image with tissue
image	
Mix map:	7 types, available when Mix
mode is active	
Persistence:	8 steps
Dynamic range:	same as B mode
Gray map:	8 types
Tint map:	Off; 8 types
Supports U/D Flip and L/R Flip	
Rotation:	same as B
CEUS Position:	On/Off
Line density:	L/M/H/UH
FOV:	On/Off
FOV size/position:	continuously adjustable
ExFov:	Off, 1-2, 1/step
Gain:	0-100, 1/step
iTouch:	On/Off, -8~8, 2/step
Image quality:	3 levels
Smooth:	0-6, 1/step
Enhance:	On/Off
Markline:	On/Off
LGC:	8 points

\*The I9 series is designed for compatibility with commercially available ultrasound contrast agents. Because the availability of these agents is subject to government regulation and approval, product features intended for use with these agents may not be commercially marketed nor made available before the contrast agent is cleared for use. Contrast related product features are enabled only on systems for delivery to an authorized country or region of use. Mindray medical systems makes no claims concerning the safety or effectiveness of contrast agents.

### Contrast Imaging QA

Support Time-Intensity Curve analysis	
Table display:	display data in table

## Performance Specifications

### Imaging (continued)

Up to 8 ROIs	
Delete all	
Delete current	
Copy ROI	
Fit curve	
Raw curve	
Motion tracking:	Reduce the effect of tissue movement
X scale:	1-5, 1/step
Auto play:	Stop, X1/10, X1/5, X1/4, X1/3, X1/2, X1, X2, X3
ROI Type:	Trace ROI, Ellipse ROI
Export:	export current data as CSV format file

### LVO

Only available in LVO exam mode	
Dedicated left ventricle contrast imaging tool	
Fetal heart:	Stomach, 4 chamber, LVOT

### Volume CEUS

Available on SD8-1s and DE11-3Ws probes	
Display:	Contrast, Tissue
Capture 3D image	
Other parameters are the same as those of 4D imaging	

### STE Imaging (Sound Touch Elastography Imaging)

Available on SC6-1s, C6-2Gs, L9-3s, L20-5s, L14-3Ws, V11-3Hs, DE11-3Ws probes	
Display Format:	V1:1, H1:1, Full
Invert:	On/Off
HQ Elasto:	On/Off
E Quality:	Pen, Gen, Res
Elas.Metric:	E, Cs, G
Scale:	30 levels
Opacity:	0-5, 1/step
Map:	3 types
ROI Width/Height:	continuously random adjustable
ROI Center Depth:	continuously adjustable
Depth:	same as B mode
iLayering:	On/Off
Filter:	0, 1
RLB View:	On/Off
M-STB Index:	On/Off
M-STB Sensibility:	0~4, 1/step
Smooth:	0~2, 1/step
Persistence:	0~2, 1/step
RLB Map:	On/Off, RLB, RLB&E, RLB&B&E
Map Position:	0%~100%, 5%/step
E bar:	Mean, Max, Min, SD
E Avg:	Off, 8 levels
Select/Bad:	On/Off
Fixed ROI:	On/Off
Save All:	On/Off
Lesion:	Off, 1-10

### Natural Touch Elastography

Available on L9-3s, L14-3Ws, L20-5s, V11-3Hs, DE11-3Ws probes	
---	--

Support strain, strain ratio and strain histogram measurement	
Unique shell analysis function	
Stress compensation technology reduces deeper tissue artifacts, obtain more uniform stress throughout whole field	
Stress indicator:	supports frame by frame stress indication
Map:	6 types
Smooth:	0-5, 1/step
Opacity:	0-5, 1/step
ROI Width/height:	continuously adjustable
Invert:	On/Off
Display Format:	V1:1, H1:1, Full
Strain mode:	0~1, 1/step
Dynamic Range:	0~5, 1/step
Map Position:	0%~100%, 5%/step
Strain Scale:	0-5, 1/step
E Sensitivity:	0-5, 1/step

### High Frame Rate STE

Available on SC6-1s, C6-2Gs, L9-3s, L20-5s, L14-3Ws, V11-3Hs, DE11-3Ws probes	
Improve the image refresh rate and provide a smoother STE image	

### STQ Imaging (Sound Touch Quantification Imaging)

Available on SC6-1s, C6-2Gs, L9-3s, L20-5s, L14-3Ws, V11-3Hs, DE11-3Ws probes	
ROI Adjustment:	adjust the ROI fixed size
Elas.Metric:	E, Cs, G
E bar:	Mean, Max, Min, SD
M-STB Index:	On/Off
M-STB Sensibility:	0-4, 1/step
Filter:	0, 1
Smooth:	0-2, 1/step
Persistence:	0-2, 1/step
Map Position:	0~100%, 5%/step
Lesion:	Off, 1~10
The square height of the elasto curve represents the average value of the elasto metric for current frame	
Scale:	0-9, 1/step
E Avg:	off, 8 levels
HQElasto:	On/Off
Select/Bad:	On/Off
Fixed ROI:	On/Off
Save All:	On/Off
Lesion:	Off, 1-10

### Auto EF

Output EDV/ESV/EF/SV/CO by Simpson method	
Activated with or without ECG	
Adjustment for the border of endocardium by single point or multi points	
Adjust Frame	
Layout:	Dual/ Single
Diastole FR	
Systole FR	
Volume curve:	On/Off

### TDI QA

Dedicated quantification tool for TDI velocity, strain, strain rate analysis	
--	--

Ellipse ROI, Standard ROI	
Up to 8 of ROI	
Delete all	
Delete current	
ROI tracking:	tracking ROI along with cardiac movement
Smooth:	1-7, 1/step
X scale:	1-5, 1/step
Std.Height:	1.5-50 mm
Std.Width:	1.5-50 mm
Std.Angle:	-89-90 degrees
Export:	export current data as CSV format file

### TT QA

Available on the SP5-1s probe under the cardiac mode	
Tissue tracking quantitative analysis	
Mandatory ECG connection before TTQA cine acquisition	
Six views for analysis:	ALAX, A4C, A2C, PSAXB, PSAXM, PSAXAP
Reload:	reload cine again for new study
Edit:	modify trace points
Start tracking	
Accept & compute:	start tracking myocardium movement when user accept trace result
Display effect:	0/1; at 1, tracking in velocity vector arrow; at 0, tracking in dots
Trace method:	3 point or manual for ALAX, A4C, A2C; manual for PSAXB, PSAXM, PSAXAP
Bull's eye:	trace result in bull's eye model
Valve's open and close time index:	MVC, MVC', AVC, AVO, MVO
Data export:	export data in CSV file
Cycle:	ECG triggered cardiac cycle recognition for analysis; adjustable
Auto play:	stop, X1/10, X1/5, X1/4, X1/3, X1/2, X1, X2, X3
Thickness:	1-30 mm, 1 mm/step; adjust trace thickness
Track point:	20-40, 1/step
Parameter:	Volume, Speed, Displacement, L Strain, L Strain R, T Strain, T Strain R, Area, R Strain, R Strain R, C Strain, C Strain R, C Rotation, C Rotation R
Smooth:	0-4, 1/step
Trace method:	3 point, manual
Tracking cycles:	1, 3
Select Cycle:	select among 3 cycles when the Tracking Cycle is set to 3
Average Cycle:	On/Off
Velocity scale:	0%-200%
Display style:	All, Endo., Myo., Epi., Curve Display
Torsion & Torsion Rate Curve	
LGC adjustment	

## Performance Specifications

### Imaging (continued)

Data Export: export current data as CSV format file

### Stress Echo

Available on the SP5-1s probe in cardiac mode  
 14 factory protocols  
 User-defined protocols  
 ECG triggered acquisition, display, selection, comparison, evaluation and archiving of multiple cardiac loops during various stages of a stress echo examination  
 Customized stages: up to 7 views per stage, and up to 12 stages per study  
 View: standard views (PLAX, SAB, PSAX, SAA, A4C, A2C, ALAX), and customized views

### Image acquisition

R-wave trigger  
 Acquire mode: Manual ROI or full screen  
 Ability to acquire frames or clips in B-mode, M-mode, Color, PW and TDI

### Image selection

Attach the images with view annotation label PLAX, SAB, PSAX, SAA, A4C, A2C, ALAX, and customized views

### Review

Automatically adjust to the number of images user defined

### Wall Motion Scoring

ASE 16 (with score 4-7), or ASE 17 (with score 4-7)  
 Graphical display of scoring (Normal, Hyperkinetic, Severely Hyperkinetic, Akinetic, Dyskinetic)

### LV volume measurement

Measurement of LV Volume in all phases of cardiac cycle

### Report

Reporting for both Wall Motion  
 Scoring and LV volume measurement

### Smart Pelvic

Enter smart pelvic in 2D or 3D/4D scanning mode  
 Set Rest, Valsalva, and Contraction frame  
 Measure automatically

### iScanHelper

Tutorial function as a guidance to show basic scanning skill with graphic of probe position, schematic of anatomy and example clinical image  
 Support ABD, SMP, URO, OB, GYN applications  
 Support broadcasting the scanning skill in multi languages

### iCompare

Allow to compare real-time ultrasound imaging to the past DICOM CT/MRI/Mammography/XRay/ Ultrasound images without external workstation  
 Helpful to easily evaluate and follow up the progression of disease, treatment effect monitoring

### DVR

Digital video recorder, a useful tool for education and memory  
 Max storage length each time: 60 min

### Cine Review and Raw Data Processing

#### Cine Review

Available in all modes  
 Frame by frame manual cine loop review or auto playback with variable speed  
 Maximum cine memory up to 63575 frames (B storage server) or 210.65 s (M storage server) (depends on the mode)  
 Maximum 4D cine memory: 5423 volumes (SD8-1s)  
 Retrospective storage (1-120s presettable) and prospective storage (1-480s pre-settable)  
 Frame compare: displays one cine in dual format and allows frame by frame compare side by side  
 Cine compare: compare cines which are saved in same imaging mode  
 Jump to first and jump to last: one keystroke go to first or last frame in the cine

#### Raw Data Processing

B-mode:  
 TGC  
 Gain  
 Dynamic range  
 Gray map  
 Tint map  
 iClear  
 L/R Flip  
 U/D Flip  
 Rotation  
 iTouch  
 LGC  
 Dual live  
 Auto Merge  
 H Scale  
 Echo Boost  
 B/iNeedle  
 Smooth  
 Zoom  
 Ref Lines  
 Dehaze  
 V1:1  
 ExtImage  
 Edge Enhance  
 M-mode:  
 Gain  
 Speed  
 Dynamic Range  
 Gray Map  
 Tint Map  
 Display format  
 Color:  
 Gain  
 Baseline  
 Smooth  
 Color map  
 Priority  
 Dual Live

Invert  
 Velocity tag  
 Glazing flow

#### PW:

Gain  
 Baseline  
 Volume  
 Angle  
 Speed  
 Dynamic range  
 Gray map  
 Tint Map  
 Display format  
 Invert  
 WF  
 T/F Res

#### Report\*

##### IMT

Intima-Media Thickness Measurement  
 Automatic detection of IMT when ROI is set  
 Support CCA, ICA, ECA, Bulb IMT  
 Near wall and far wall detection  
 Angle selectable

##### Smart OB™

Auto measurement for OB, a special too for easy OB scan and greatly reduce time and increase productivity  
 Support BPD, HC, OFD, FL, AC  
 Better get GA before start auto AC  
 Measurement result can be modified by user

##### Report

Specific report template by application  
 Editable value in report  
 Images selectable  
 Anatomy information for vascular and OB report  
 Editing though iReport  
 User-defined report template  
 Selecting report modules  
 Adding/removing measurement items from the report  
 Changing report layout  
 Load/save comment  
 Viewing history reports  
 Preview and printing reports  
 Able to Export as PDF/RTF file  
 V-Mapping  
 Add anatomical graphics for illustration

##### Mini report

Quickly displaying Mini report in the thumbnail area of the main screen  
 Including both general measurement and application measurement results  
 Support deleting measurement results

## Performance Specifications

### Report\* (continue)

#### iWorks

Auto workflow protocol	
Templates are user configurable	
Functions:	pause, stop, replace, repeat, skip, insert single step, return and continue, steps in thumbnail, iNSert™ another template
iWorks setup mode:	B; B/B (Dual Live); Dual B/B; Color; B/Color (Dual Live); Power; B/Power (Dual Live); B + PW; Color + PW; Power + PW; B + CW; Color + CW; B+M; B+TVI; TVI+TVD; iScape.
iWorks setup annotation:	support up to 2 annotations, location and font size are configurable.
iWorks setup bodymark:	select existing library, and transducer indicator is pre-settable
iWorks setup measurement:	select existing measurement library
Template import and export are available	
Support create user-defined iWorks protocol	

#### UltraView™

Components:	
DICOM Basic	
DICOM Query/Retrieve	
DICOM OB/GYN SR	
DICOM Cardiac SR	
DICOM Vascular SR	
DICOM Breast SR	
DICOM Urology SR	
DICOM Pediatric SR	
DICOM Small Parts SR	
DICOM Abdomen SR	
TDI QA	
Contrast Imaging QA	
Tissue Tracking QA	
Stress Echo	
SCV (Slice Contrast View)	
iLive	
iPage+	
IVF	
Ultrasound Fusion Imaging	
Auto EF	
Smart Planes CNS	
PC-based off-line software	

\*Not all measurements are listed in this part; For more detailed information please refer to User Manual

### Exam Storage and Management

#### Exam Storage

1TB hard drive
128G SSD (Solid State Drive)
Direct digital storage of single frame and cine 2D,

color and Doppler

#### Exam Management

iStation™ workstation dedicated for patient exam management
Patient exam query/retrieve
Support review of current and past exam
New exam, Activate exam, Continue exam functions, End exam are available
Support measurements and calculations on archived exam and images
Export images in BMP/JPG/TIFF/DCM/AVI/MP4 format Support backup/send to USB devices, DVD-RW, CD-R, DVD+R, DVD-R, DVD+RW media

#### Connectivity

##### Ethernet Network Connection

Cable connection
Wireless connection: built-in wireless adapter

##### DICOM 3.0

DICOM basic	
Verify (SCU, SCP)	
Print	
Store	
Storage Commitment	
Media Exchange	
DICOM Worklist	
DICOM Query/Retrieve	
DICOM Modality Performed Procedure Step	
MPPS	
DICOM OB/GYN structure report	
DICOM Cardiac structure report	
DICOM Vascular structure report	
DICOM Breast Report	
DICOM Abdominal structure report	
DICOM Small Parts structure report	

##### iStorage (included in UltraAssist)

Direct network storage tool between ultrasound system and personal computer
---

##### MedSight

An interactive app that lets you transfer clinical images straight from Mindray Ultrasound system to a smart device, such as mobile phone or tablet PC	
Needs to be installed on mobile terminal	
Transfer images or clips from system to mobile terminal through Wi-Fi	
Support both iOS (7.0 and above) and Android (4.0 and above) system	
For iOS powered smart device:	DICOM is mandatory
For Android powered smart device:	DICOM not necessary

##### MedTouch

Connect Ultrasound machine to smart devices based on Android and iOS system, such as tablet PC or mobile phone. Remote control of Ultrasound machine, review of patient information, and tutorial software iScanHelper study on smart devices
---

Support Android and iOS powered smart devices
Android 4.0 and above
iOS 7.0 and above
DICOM not necessary

#### Transducers

##### Curved Array

SC6-1s	
Application:	Obstetrics, Gynecology, Abdomen, Small Parts, Musculoskeletal, Vascular, Urology, Nerve
Bandwidth:	1.2-6.0 MHz
Number of Elements:	192
FOV (max):	60°
Extended FOV:	72°
Convex Radius:	60 mm
Depth:	4-40 cm
Physical Footprint:	65.1 mm × 16.4 mm
Footprint:	64.9 mm × 16.2 mm
B-mode	
Frequencies:	1.2~3.8, 1.7~5.2, 2.0~6.0 MHz
Harmonic	
Frequencies:	4.0, 5.0, 6.0 MHz
Color Frequencies:	2.0, 2.5, 3.0, 3.3 (HR Flow) MHz
PW Frequencies:	2.0, 2.5, 3.0 MHz
Biopsy Guide:	NGB-022, multi angle, reusable; C11-3s
Application:	Abdomen, Vascular, Cardiac, Small Parts, Pediatric, Cephalic
Bandwidth:	2.6-12.8 MHz
Number of Elements:	128
FOV (max):	101°
Extended FOV:	113°
Convex Radius:	15 mm
Depth:	1.5-35 cm
Physical Footprint:	32.8 mm × 25 mm
Footprint:	27.4 mm × 8.4 mm
B-mode	
Frequencies:	2.6~6.5, 3.2~7.9, 4.7~12.8 MHz
Harmonic	
Frequencies:	6.0, 7.0, 8.0 MHz
Color Frequencies:	4.4, 5.0, 5.7, 5.0 (HR Flow) MHz
PW Frequencies:	4.4, 5.0, 5.7 MHz
Biopsy Guide:	NGB-018, multi angle, reusable

##### C6-2Gs

Application:	Obstetrics, Gynecology, Abdomen, Urology
Bandwidth:	1.2-6.0 MHz
Number of Elements:	128
FOV (max):	94°
Extended FOV:	106°
Convex Radius:	20 mm
Depth:	4-40 cm
Physical	

## Performance Specifications

### Transducers (continued)

Footprint:	37.6 mm × 19 mm
Footprint:	31.5 mm × 11.2 mm
B-mode	
Frequencies:	1.2~3.8, 1.7~5.2, 2.0~6.0 MHz
Harmonic	
Frequencies:	3.2, 4.0, 5.0 MHz
Color Frequencies:	2.0, 2.5, 3.0, 3.5 (HR Flow) MHz
PW Frequencies:	2.0, 2.5, 3.0 MHz
Biopsy Guide:	NGB-024, multi angle, reusable;

### V11-3Hs

Application:	Obstetrics, Gynecology, Urology
Bandwidth:	3.0-11.0 MHz
Number of Elements:	192
FOV (max):	170°
Extended FOV:	210°
Convex Radius:	11 mm
Depth:	1.5-28 cm
Physical Footprint:	24.9 mm × 21.8 mm
Footprint:	24 mm × 9 mm
B-mode	
Frequencies:	3.0~7.0, 4.0~9.0, 5.0~11.0 MHz
Harmonic	
Frequencies:	8.0, 9.0, 10.0 MHz
Color Frequencies:	4.4, 5.0, 5.7, 5.5 (HR Flow) MHz
PW Frequencies:	4.5, 5.0, 5.5 MHz
Biopsy Guide:	NGB-025, single angle, reusable

### SC8-2s

Application:	Obstetrics, Gynecology, Abdomen, Urology, Vascular
Bandwidth:	1.8-8.2 MHz
Number of Elements:	192
FOV (max):	76°
Extended FOV:	88°
Convex Radius:	40 mm
Depth:	4-40 cm
Physical Footprint:	26.3 mm × 66.9 mm
Footprint:	15 mm × 52 mm
B-mode	
Frequencies:	1.8~5.4, 2.3~6.8, 2.8~8.2 MHz
Harmonic	
Frequencies:	4.0, 5.5, 6.0 MHz
Color Frequencies:	3.0, 3.5, 4.0, 3.5 (HR Flow) MHz
PW Frequencies:	3.0, 3.5, 3.8 MHz
Biopsy Guide:	NGB-029, multi angle, reusable

### Volume Curved Array

#### SD8-1s

Application:	Obstetrics, Gynecology, Abdomen
Bandwidth:	1.8-8.2 MHz
Number of Elements:	192
FOV (max):	66°
Extended FOV:	91°
Volume Sweep	
Angle (max):	85°

Convex Radius:	45 mm
Depth:	4-40 cm
Physical Footprint:	75.7 mm × 52.6 mm
B-mode	
Frequencies:	1.8-5.4, 2.3-6.8, 2.8-8.2 MHz
Harmonic	
Frequencies:	4.0, 5.0, 5.5, 6.0 MHz
Color Frequencies:	3.0, 3.5, 4.0, 3.5 (HR Flow) MHz
PW Frequencies:	3.0, 3.5, 3.8 MHz
Biopsy Guide:	NGB-039, multi angle, reusable

### DE11-3Ws

Application:	Obstetrics, Gynecology, Urology
Bandwidth:	2.0-9.0 MHz
Number of Elements:	192
FOV (max):	162°
Extended FOV:	187°
Volume Sweep	
Angle (max):	120°
Convex Radius:	10 mm
Depth:	4-40 cm
Physical	
Footprint:	24 mm × 24 mm
B-mode	
Frequencies:	2.0~6.0, 2.8~8.2, 3.0~9.0 MHz
Harmonic	
Frequencies:	4.0, 5.0, 6.0 MHz
Color Frequencies:	4.4, 5.0, 5.7, 5.0 (HR Flow) MHz
PW Frequencies:	4.4, 5.0, 5.7 MHz
Biopsy Guide:	NGB-047, single angle, reusable

### Linear

#### L9-3s

Application:	Abdomen, Pediatric, Small Parts, Musculo-skeletal, Vascular, Nerve, Obstetrics
Bandwidth:	2.5-9.0 MHz
Number of Elements:	192
Field of View (max):	43.7 mm
Extended FOV:	60° (OB1, NT); 40° (Others)
Steered Angle:	B: 0°, +/-6°, +/-12°; C/PW: -30°-30°
Depth:	1.5-35 cm
Physical Footprint:	62 mm × 22 mm
Footprint:	48 mm × 11 mm
B-mode	
Frequencies:	2.5~7.0, 3.4~8.2, 3.6~9.0 MHz
Harmonic	
Frequencies:	5.0, 6.0, 7.0 MHz
Color Frequencies:	3.0, 3.6, 5.0, 4.0 (HR Flow) MHz
PW Frequencies:	3.0, 3.6, 5.0 MHz
Biopsy Guide:	NGB-034, multi angle, reusable

#### L20-5s

Application:	Abdomen, Small Parts, Musculo-skeletal, Vascular, Nerve
Bandwidth:	6-23 MHz
Number of	

Elements:	192
Field of View (max):	28.5 mm
Extended FOV:	20°
Steered Angle:	B: 0, +/-6°, +/-12°; C/PW: -20°-20°

Depth:	1.5-29 cm
Physical Footprint:	42.23 mm × 22.10 mm
Footprint:	31.5 mm × 4.5 mm
B-mode	
Frequencies:	6.0~13.0, 9.0~16.6, 12.5~23.0 MHz

Harmonic	
Frequencies:	14.0, 16.0, 18.0 MHz
Color Frequencies:	9.0, 11.0, 13.0, 13.0 (HR Flow) MHz
PW Frequencies:	8.3, 10.0, 12.5 MHz
Biopsy Guide:	not available

### L14-3Ws

Application:	Small Parts, Musculoskeletal, Vascular, Abdomen, Pediatric, Thoracic/Pleural
Bandwidth:	3.0-14.0 MHz

Number of Elements:	256
Field of View (max):	50.8 mm
Extended FOV:	20°
Steered Angle:	B: 0, +/-6°, +/-12°; C/PW: -30°-30°

Depth:	1.5-35 cm
Physical Footprint:	66.8 mm × 25.5 mm
Footprint:	55.5 mm × 8.2 mm
B-mode	
Frequencies:	3.0~9.0, 5.0~12.0, 6.0~14.0 MHz

Harmonic	
Frequencies:	10.0, 11.0, 12.0, 13.0 MHz
Color Frequencies:	5.0, 6.2, 7.2, 8.3 (HR flow) MHz
PW Frequencies:	5.0, 6.2, 7.2 MHz
Biopsy Guide:	NGB-054, multi angle/depth, reusable

### Phased Array

#### P7-3Ts

Application:	Cardiac
Bandwidth:	2.3-7.2 MHz
Number of Elements:	64
Field of View (max):	90°
Extended FOV:	90°
Depth:	2.0-38.0 cm
Physical Footprint:	14 mm × 12 mm
Footprint:	12.2 mm × 12.2 mm
B-mode	
Frequencies:	2.3~5.4, 2.8~6.4, 3.3~7.2 MHz
Harmonic	
Frequencies:	5.0, 6.0, 7.0 MHz
Color Frequencies:	2.7, 3.3, 4.0, 4.0 (HR Flow) MHz; TDI: 5.0, 6.2 MHz
PW Frequencies:	2.7, 3.3, 4.0 MHz
Biopsy Guide:	not available

## Performance Specifications

### Transducers (continued)

#### SP5-1s

Application:	Abdomen, Cardiac, Vascular, Cephalic, Thoracic/Pleural
Bandwidth:	1.5-4.5 MHz
Number of Elements:	80
Field of View (max):	90°
Extended FOV:	90°
Depth:	2-38 cm
Physical Footprint:	38.2 mm × 30.5 mm
Footprint:	23.4 mm × 15.2 mm
B-mode Frequencies:	1.5~2.5, 2.5~3.5, 3.5~4.5 MHz
Harmonic Frequencies:	3.0, 3.5, 3.5, 4.0, 4.0 MHz
Color Frequencies:	2.0, 2.3, 2.5, 2.5(HR Flow) MHz; TDI: 3.0, 3.8 MHz
PW Frequencies:	2.0, 2.3, 2.5 MHz; TDI: 2.5, 4.0 MHz
CW Frequencies:	2.0 MHz
Biopsy Guide:	NGB-011, multi angle, reusable

#### Pencil

##### CW5s

Application:	Vascular, Cephalic, Pediatric
Number of Elements:	2
Biopsy Guide:	not available

##### CW2s

Application:	Cardiac, Cephalic, Pediatric
Number of Elements:	2
Biopsy Guide:	not available

### Peripheral Devices and Accessories (Option)

#### Black/White Video Printer

Digital	MITSUBISHI P95DW-N
Analog	SONY UP-X898MD

#### Color Digital Printer

SONY UP-D25MD

#### Graph/Text Printer

HP OFFICEJET PRO 8100

#### Gel Warmer

Support gel warming with 3 angle position:	15, 55, 90 degrees
--	--------------------

Easily be disassembled off system for cleaning

Temperature with 4 levels:	off/34° C/37° C/40° C, with deviation of ±1° C
----------------------------	--

Light indicator for temperature protecting	
Dimension:	82(D)*78(W)*119(H) mm
Weight:	approx. 240 g (net)

Continuous operation time:	>12 h
----------------------------	-------

#### Footswitch

USB port:	FS-81-SP-2 (single pedal), 971-SWNOM (2/3-pedal)
-----------	--

Support User-definable functions (Freeze, Save, Print)

#### ECG

6-pin, AHA/IEC, for 3-lead wires	
ECG wave display:	on/off
ECG source:	Lead/External
Position:	0~100%, 5%/step
Trig mode:	off/single/dual/timer
Gain:	0-30, 1/step
Sweep speed:	6 steps
Invert:	on/off

#### PCG (not for sale in EU countries)

PCG wave display:	on/off
Gain:	0-30, 1/step
Speed:	6 steps
Smooth:	1-4, 1/step

#### Barcode Reader

SYMBOL LS2208 (1D)	
SYMBOL DS4308 (2D)	

#### Built-in Wi-Fi 5 Wireless Adapter

Encryption:	WPA, WPA2
Max transfer speed:	300 Mbps
Protocols:	IEEE 802.11 ac/a/b/g/n
Frequency:	2.4G/5G

#### iVocal Microphone

SAMSON XPD2

#### Built-in Battery

Replaceable and rechargeable lithium battery	
Full battery lasts for more than 24h in standby mode	
Battery capacity indicators without powering on the system	
Battery fully-recharged time:	less than 6 h
Continuous scanning time:	more than 120 mins for 2H battery

### Probe Adapter

PCM-ES01:	transforming E socket to S socket, only for SC6-1E, SD8-1E, P7-3TE, and SP5-1E
PCM-US01:	transforming U socket to S socket, only for P7-3TU and SP5-1U

### System Inputs and Outputs

#### Audio input/Output

Microphone/Audio:	1 port
-------------------	--------

#### Video Output

S-Video out:	1 port, PAL/NTSC
HDMI:	1 Port
VGA out:	1 port

#### Physio input

Support ECG/PCG signal	
ECG:	1 port
PCG:	1 port

#### Other input/Output

USB:	6 ports (5 USB 3.0 and 1 Type-C)
Ethernet:	1 port

### Safety and Conformance

#### Quality Standards

ISO 9001	
ISO 13485	

#### Design Standards

EN 60601-1 and IEC 60601-1	
EN 60601-1-2 and IEC 60601-1-2	
EN 60601-1-6 and IEC 60601-1-6	
EN 60601-2-37 and IEC60601-2-37	
EN 62304 and IEC 62304	
EN 62366 and IEC 62366	
EN ISO 17664 and ISO 17664	

#### 12.3 CE Declaration

The ultrasound system is fully in conformance with the Council Directive 93/42/EEC Concerning Medical Devices. The number adjacent to the CE marking (0123) is the code of the EU-notified body that certified meeting the requirements of Annex II excluding (4). of the Directive.

#### NOTICE:

Not all features or specifications described in this document may be available in all transducers and/or modes. Mindray reserves the right to make changes in specifications and features shown herein, or discontinue the product at any time without notice or obligation.

### Mindray North America

800 MacArthur Boulevard, Mahwah, NJ 07430

Tel: 800.288.2121 Support: 877.913.9663 Fax: 800.926.4275 [www.mindray.com](http://www.mindray.com)