

# DP-50 Vet

Hand-carried Diagnostic Ultrasound System

## Exceptional ergonomics by design

- Sleek, streamlined, compact shape
- Lightweight system for enhanced mobility
- High resolution, wide-angle 15" LCD with tilt functionality for better viewing

## Diagnostic confidence with advanced imaging technologies

- iTouch™ Auto Optimization: One touch image optimization, quickly optimizes the uniformity and dynamic range, very convenient and efficient
- iClear™ Speckle Suppression: Scale the image size to intelligently detect structures of different sizes, improve accuracy
- Phase Shift Harmonic Imaging: Enhance contrast resolution and improve image quality by separating basal frequency and harmonic content

## Efficient and intelligent workflow

- Quick Save: One push image transfer to local or directly to USB
- Onboard reporting: Standard templates for various applications
- DICOM 3.0 Compliant: Structured reporting, Worklist, DICOM print, Storage



# DP-50 Vet

Hand-carried Diagnostic Ultrasound System



## Exceptional service and support

Mindray is a global medical device business with three established segments: Patient Monitoring and Life Support Products, *In-Vitro* Diagnostic Products and Medical Imaging Systems. With more than 6,000 employees and with R&D centers on three continents, Mindray offers the medical community a range of innovative solutions designed to contain costs while ensuring the highest quality of patient care.

Mindray's clinical education specialists provide your staff with comprehensive, on-site training and in-service on demand 24-hours a day, 7-days a week.

Our dedicated field service and in-house technical support organization are one of the largest direct service teams in the industry. Individually and collectively they are committed to enhancing patient care while protecting your investment.

Learn how we can help you maximize your resources with some of the most advanced and affordable devices in the industry. For more information about Mindray products and services, visit [www.na.mindray.com](http://www.na.mindray.com).



Mindray DS USA, Inc.  
U.S. headquarters in Mahwah, NJ

Mindray DS USA, Inc.  
800 MacArthur Blvd., Mahwah, NJ 07430  
Tel: 1.800.288.2121 Tel: 201.995.8000 Fax: 1.800.926.4275 [www.na.mindray.com](http://www.na.mindray.com)

Mindray™ is a trademark or registered trademark of Shenzhen Mindray Bio-Medical Electronics Co., Ltd.  
© 2012 Mindray DS USA, Inc. Subject to change. P/N: 0002-08-45003 Rev A

**mindray**  
NORTH AMERICA