



iPM12 Vet iMEC8 Vet

Safety & Efficiency
for Veterinary Monitoring



mindray




iPM12 Vet

Safety & Efficiency for Veterinary Monitoring

- 12.1" high resolution touchscreen display with up to 8 waveforms
- Lightweight, ergonomic design with incorporated carry handle
- Standard features include 3 or 5 lead ECG, 2 IBP, SpO₂, NIBP, respiration, temperature, integrated recorder
- Additional module parameter options include:
 - Sidestream or Mainstream CO₂
 - Multi-gas
- Extensive data storage capability for trend data, alarms, events, NIBP measurements and up to 48 hours of full disclosure
- Integrated 3-trace recorder

360-degree visible alarm indicator

Plug-and-play modules for several different combinations of advanced parameters: IBP, EtCO₂, and Multi-gas



Customizable quick keys provide quick, convenient access to frequently used functions

Rotary knob and hard keys for convenient operation




iMEC8 Vet

Safety & Efficiency for Veterinary Monitoring

- 8.4" high-resolution LED display with up to 8 waveforms
- Compact portable design
- Trust for veterinary use, parameters including 3 or 5 lead ECG, IBP, SpO₂, NIBP, respiration, temperature, integrated recorder
- Additional parameter options include:
 - Sidestream CO₂
 - Mainstream CO₂
- Up to 4 hours continuous monitoring with large capacity Li-ion battery
- 360-degree visible alarm indicator
- 120 hours of trend data and up to 48 hours of full disclosure

Built-in handle for ease-of-use during transport

Li-ion battery allows for up to 4 hours of continuous monitoring



Rotary knob and hard keys for convenient operation

Integrated 3-trace thermal recorder provides waveform and data reports at point-of-care



- Mindray Global Offices
- ◉ Mindray R&D / Mfg Operations


Mindray is a global medical device manufacturer with three established business segments: Patient Monitoring and Life Support Products, Medical Imaging Systems, and *In-Vitro* Diagnostic Products. With more than 7,000 employees and R&D centers on three continents, we offer the medical community a wide range of innovative solutions designed to ensure the highest quality of care, while containing costs.

Learn how we can help you maximize your resources, with some of the most advanced and cost effective solutions in the industry.

For more information about Mindray products and services, visit www.mindray.com.



Mindray North America Headquarters
Mahwah, New Jersey



Mindray DS USA, Inc.
800 MacArthur Blvd., Mahwah, NJ USA 07430
Tel: 800.288.2121 or 201.995.8000 Fax: 800.926.4275 www.mindray.com

Mindray® is a registered trademark of Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All brands and product names are trademarks of their respective owners. ©2018 Mindray DS USA, Inc. Subject to change. 02/18 P/N: 0002-08-8027 Rev C

mindray
NORTH AMERICA